

Omnifit® Series II Cup Inserts

Catalog Number	Description	Inner Diameter	GAP Outer Diameter	Trial Catalog Number
2041-XX46	10°	22, 26, 28mm	48mm	2205-XX46
2041-XX50	10°	22, 26, 28, 32mm	52mm	2205-XX50
2041-XX54	10°	22, 26, 28, 32mm	56mm	2205-XX54
2041-XX58	10°	22, 26, 28, 32mm	60mm	2205-XX58
2041-XX62	10°	22, 26, 28, 32mm	64mm	2205-XX62
2041-XX66	10°	22, 26, 28, 32mm	68mm	2205-XX66
2041-XX70	10°	22, 26, 28, 32mm	72mm	2205-XX70
2042-XX46	20°	22, 26, 28mm	52mm	2206-XX46
2042-XX50	20°	22, 26, 28, 32mm	52mm	2206-XX50
2042-XX54	20°	22, 26, 28, 32mm	56mm	2206-XX54
2042-XX58	20°	22, 26, 28, 32mm	60mm	2206-XX58
2042-XX62	20°	22, 26, 28, 32mm	64mm	2206-XX62
2042-XX66	20°	22, 26, 28, 32mm	68mm	2206-XX66
2042-XX70	20°	22, 26, 28, 32mm	72mm	2206-XX70
2043-XX46	0°	22, 26, 28mm	48mm	2206-XX46
2043-XX50	0°	22, 26, 28, 32mm	52mm	2213-XX50
2043-XX54	0°	22, 26, 28, 32mm	56mm	2213-XX54
2043-XX58	0°	22, 26, 28, 32mm	60mm	2213-XX58
2043-XX62	0°	22, 26, 28, 32mm	64mm	2213-XX62
2043-XX66	0°	22, 26, 28, 32mm	68mm	2213-XX66
2043-XX70	0°	22, 26, 28, 32mm	72mm	2213-XX70

Restoration™ GAP II Acetabular Shell

Catalog Number	Description
2082-0048L,R	48mm OD
2082-0052L,R	52mm OD
2082-0056L,R	56mm OD
2082-0060L,R	60mm OD
2082-0064L,R	64mm OD
2082-0068L,R	68mm OD
2082-0072L,R	72mm OD

Restoration™ GAP Plate Screws

Catalog Number	Description
2080-00XX	6.5mm diameter 15-60mm length (in 5mm increments)

Instrumentation

Instrument	Catalog Number
Mini Bone Plate Bender	2107-5004
Inferior Hook Crimper	2107-5010
Superior Plate Cutter	2107-5020

Important:

- Use only Restoration™ GAP Plate Screws (Series 2080) for all screw holes
- Avoid repeated bending and sharp notch for inferior hook and superior plates
- The position of the obturator artery should be identified and avoided prior to the placement of the inferior hook
U.S. Patent 5,702,477

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

The marks bearing the symbol "TM" are trademarks of Stryker.
The marks bearing the symbol "R" are registered trademarks of Stryker.

Literature Number: LRST-2
GC/GS 900 01/05

Copyright © 2004 Stryker
Printed in USA

stryker®

325 Corporate Drive
Mahwah, NJ 07430
t: 201 831 5000

www.stryker.com