



**The information contained in this document
is intended for healthcare professionals only.**



EXETER™
total hip system
V40™

stryker® **Howmedica** OSTEONICS

Size Options



Femoral Head Offset (mm)

		-4 Head	+0 Head	+4 Head	+8 Head
Stem	30 (AP1)	27	30	33	36
	33 (AP2)	30	33	36	39
	35 (Small CDH)	32	35	38	41
	37.5 (#0-4)	34	37.5	40.5	43.5
	44 (#0-5)	41	44	47	50
	50 (#1-5)	47	50	53	56

Leg Length Increase (mm)*

		-4 Head	+0 Head	+4 Head	+8 Head
Stem	30 (AP1)	-2	-	+2	+4
	33 (AP2)	-2	-	+2	+4
	35 (Small CDH)	-2	-	+2	+4
	37.5 (#0-4)	-2	-	+2	+4
	44 (#0-5)	-2	-	+2	+4
	50 (#1-5)	-2	-	+2	+4

* Assume Middle Reference Hole

Stryker® is a registered trademark of Stryker Corporation.
Exeter™ and V40™ are trademarks of Howmedica Osteonics Corp.
©Howmedica Osteonics Corp. 2002 Lit. No. **LEV40WC** 1.6M 05/02 MC/ITP Printed in U.S.A..

stryker®
Howmedica
OSTEONICS

59 Route 17 South
Allendale, NJ 07401
howost.com