

**The information contained in this document
is intended for healthcare professionals only.**

Trident[®]
ALL POLY ACETABULAR SYSTEM



Trident[®]
ALL POLY ACETABULAR SYSTEM

The first all poly cup available in the market featuring highly crosslinked polyethylene.



CROSSFIRE[®]




stryker[®]
Howmedica
OSTEONICS

ELEVATED RIM
Designed to Improve
Stability

**CROSSFIRE® HIGHLY
CROSSLINKED POLY**
Reduces gravimetric wear

X-RAY WIRE
Easy evaluation of cup location

CEMENT PRESSURIZING FLANGE
Uniform cement placement



PMMA SPACERS
Allows a minimum
2mm cement mantle

X-RAY BEAD
Easy evaluation of cup location

**STABILIZING RADIAL
PROJECTIONS**
Designed to Improve Stability

Crossfire®		Min. Poly Thickness (mm)	Trials
Catalog Number	Description		
65C-2240	22 x 40	6.4	1030-2240
65C-2242	22 x 42	7.4	1030-2242
65C-2244	22 x 44	7.9	1030-2244
65C-2642	26 x 42	5.5	1030-2642
65C-2644	26 x 44	6.5	1030-2644
65C-2646	26 x 46	7.5	1030-2646
65C-2648	26 x 48	8.0	1030-2648
65C-2844	28 x 44	5.5	1030-2844
65C-2846	28 x 46	6.5	1030-2846
65C-2848	28 x 48	7.5	1030-2848
65C-2850	28 x 50	8.0	1030-2850
65C-2852	28 x 52	8.5	1030-2852
65C-2854	28 x 54	9.5	1030-2854
65C-2856	28 x 56	10.5	1030-2856
65C-2858	28 x 58	11.5	1030-2858
65C-2860	28 x 60	12.5	1030-2860
65C-2863	28 x 63	14.0	1030-2863
65C-2866	28 x 66	15.5	1030-2866
65C-3248	32 x 48	5.5	1030-3248
65C-3250	32 x 50	6.5	1030-3250
65C-3252	32 x 52	7.5	1030-3252
65C-3254	32 x 54	8.0	1030-3254
65C-3256	32 x 56	8.5	1030-3256
65C-3258	32 x 58	9.5	1030-3258
65C-3260	32 x 60	10.5	1030-3260
65C-3263	32 x 63	12.0	1030-3263
65C-3266	32 x 66	13.5	1030-3266

Trident® All Poly Cup Acetabular Instrumentation

Cutting Edge™ Acetabular Spherical	
Reamers 36-66mm (1mm increments)	2102-04xx
Cutting Edge™ Acetabular Reamer Handle	2102-0410
Trident® All Poly Cup Introducer Handle	1030-1500
Cup Introducer Plate (22, 26, 28 & 32)	1030-16xx
Cement Pressurizer Handle	1030-1700
Window Trial Handle	2101-0200
Trident® All Poly Cup Pusher	1030-1600
Cutting Edge™ Spherical Cup Window	
Trials 40-66 (2mm increments)	2208-10xx
Trident® All Poly Cup Trials (ID, OD)	1030-xxxx
Trident® All Poly Acetabular Drill	1030-1201
Trident® All Poly Flexible Drill Shaft	1030-1200
Trident® All Poly Drill Guide	1030-1250
Trident® All Poly Cup Tray	1030-4000

- **Trident® All Poly Surgical Technique** **LTAP/ST**
- **Trident® All Poly Templates** **LTEM58B 1-2**
- **Surgical Simplex® P Bone Cement**
 - **Full Dose (10 pack)** **6191-1-010**
 - **Half Dose (10 pack)** **6188-1-010**

* Additional 22mm and 26mm sizes are now available upon special request.

trident-hip.com

Stryker® is a registered trademark of Stryker Corporation.

Trident®, Crossfire® and Simplex® P are registered trademarks of Howmedica Osteonics Corp.

Cutting Edge™ is a trademark of Howmedica Osteonics Corp.

©Howmedica Osteonics 2002 Lit. No. LTAP/SS 2M 11/02 MC/GS Rev.2 Printed in U.S.A.

stryker®
Howmedica
OSTEONICS
howost.com